


## Product Datasheet

### Mouse Abca7 ELISA Kit (orb565501)

<b>Description</b>	Mouse Abca7 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	Abca7
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q91V24</a>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Abca7(ATP-  
binding  
cassette  
sub-fam...