


Product Datasheet

Mouse C5 ELISA Kit (orb565566)

Description	Mouse C5 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	C5
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P06684
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
C5(Complement
C5) ELISA Kit