


Product Datasheet

Mouse CYP19A1 ELISA Kit (orb565585)

Description	Mouse CYP19A1 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Cyp19a1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P28649
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Cyp19a1(Cytochrome
P450 19A1) ELIS...