


## Product Datasheet

### Mouse ANXA5 ELISA Kit (orb565696)

<b>Description</b>	Mouse ANXA5 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	ANXA5
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P48036</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
ANXA5(Annexin  
A5) ELISA Kit