


Product Datasheet

Mouse CFD ELISA Kit (orb565760)

Description	Mouse CFD ELISA Kit
Reactivity	Mouse
Range	46.875-3000pg/ml
Target	CFD
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P03953
Sensitivity	28.125pg/ml
Expiration Date	Please enquire.



Mouse
CFD(Complement
Factor D) ELISA
Kit