


Product Datasheet

Mouse CHGA ELISA Kit (orb565765)

| | |
|------------------------|-----------------------|
| Description | Mouse CHGA ELISA Kit |
| Reactivity | Mouse |
| Range | 78.125-5000pg/ml |
| Target | CHGA |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P26339 |
| Sensitivity | 46.875pg/ml |
| Expiration Date | Please enquire. |



Mouse
CHGA(Chromogranin
A) ELISA Kit