


Product Datasheet

Mouse CLCF1 ELISA Kit (orb565772)

Description	Mouse CLCF1 ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	CLCF1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9QZM3
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
CLCF1(Cardiotrophin
Like Cytokine ...