


## Product Datasheet

### Mouse CORT ELISA Kit (orb565785)

<b>Description</b>	Mouse CORT ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	2.81-180ng/ml
<b>Target</b>	CORT
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sensitivity</b>	1.688ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
CORT(Corticosterone)  
ELISA Kit