


Product Datasheet

Mouse EPCAM ELISA Kit (orb565833)

| | |
|------------------------|------------------------|
| Description | Mouse EPCAM ELISA Kit |
| Reactivity | Mouse |
| Range | 0.078-5ng/ml |
| Target | EPCAM |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q99JW5 |
| Sensitivity | 0.047ng/ml |
| Expiration Date | Please enquire. |



Mouse
EPCAM(Epithelial
Cell Adhesion
Mol...