


## Product Datasheet

### Mouse FIL1d ELISA Kit (orb565850)

<b>Description</b>	Mouse FIL1d ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	FIL1d
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q9QYY1</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
FIL1d(Interleukin  
1 Delta) ELISA  
K...