


Product Datasheet

Mouse FSH ELISA Kit (orb565857)

Description	Mouse FSH ELISA Kit
Reactivity	Mouse
Range	2.344-150mIU/ml
Target	FSH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P35378
Sensitivity	1.406mIU/ml
Expiration Date	Please enquire.



Mouse FSH(Follicle-stimulating hormone) ...