


## Product Datasheet

### Mouse GLP-1 ELISA Kit (orb565892)

<b>Description</b>	Mouse GLP-1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	GLP-1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P55095</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
GLP-  
1(Glucagon  
Like  
Peptide 1)  
ELI...