


## Product Datasheet

### Mouse IL-1R1 ELISA Kit (orb565961)

<b>Description</b>	Mouse IL-1R1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	93.75-6000pg/ml
<b>Target</b>	IL-1R1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P13504</b>
<b>Sensitivity</b>	56.25pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse IL-1R1(Interleukin 1 Receptor Type...