


Product Datasheet

Mouse I-TAC ELISA Kit (orb565977)

| | |
|------------------------|------------------------|
| Description | Mouse I-TAC ELISA Kit |
| Reactivity | Mouse |
| Range | 62.5-4000pg/ml |
| Target | I-TAC |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q9JHH5 |
| Sensitivity | 37.5pg/ml |
| Expiration Date | Please enquire. |



Mouse I-TAC(Interferon Inducible T-Cell ...