


Product Datasheet

Mouse ITLN1 ELISA Kit (orb565978)

Description	Mouse ITLN1 ELISA Kit
Reactivity	Mouse
Range	0.313-20ng/ml
Target	ITLN1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	O88310
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Mouse
ITLN1(Intelectin
1/Omentin)
ELISA ...