


Product Datasheet

Mouse LTF ELISA Kit (orb565989)

Description	Mouse LTF ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/ml
Target	LTF/LF
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P08071
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Mouse
LTF/LF(Lactoferrin)
ELISA Kit