


## Product Datasheet

### Mouse MIP-3alpha ELISA Kit (orb565999)

<b>Description</b>	Mouse MIP-3alpha(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	7.813-500pg/ml
<b>Target</b>	MIP-3 $\alpha$
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<b>O89093</b>
<b>Sensitivity</b>	4.688pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse MIP-3 $\alpha$ (Macrophage Inflammatory Protein 3) Pr...