


Product Datasheet

Mouse MMP-17 ELISA Kit (orb566003)

Description	Mouse MMP-17 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	MMP-17
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9R0S3
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse MMP-17(Matrix Metalloproteinase 17...