


Product Datasheet

Mouse MMP-24 ELISA Kit (orb566004)

Description	Mouse MMP-24 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	MMP-24
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9R0S2
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse MMP-24(Matrix Metalloproteinase 24...