


Product Datasheet

Mouse MT ELISA Kit (orb566007)

Description	Mouse MT ELISA Kit
Reactivity	Mouse
Range	7.813-500pg/ml
Target	MT
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Mouse
MT(Melatonin)
ELISA Kit