


## Product Datasheet

### Mouse PAR2 ELISA Kit (orb566044)

<b>Description</b>	Mouse PAR2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	PAR2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P55086</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
PAR2(Protease  
Activated  
Receptor 2...