


Product Datasheet

Mouse Pro-ANP ELISA Kit (orb566074)

Description	Mouse Pro-ANP ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Pro-ANP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse Pro-ANP(Pro Atrial Natriuretic Pep...