


Product Datasheet

Mouse PTGS1 ELISA Kit (orb566079)

Description	Mouse PTGS1 ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	PTGS1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P22437
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
PTGS1(Prostaglandin
Endoperoxide S...