


Product Datasheet

Acetyl CoA Carboxylase (phospho-Ser79) antibody (orb5661)

Description	Acetyl CoA Carboxylase (phospho-Ser79) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Acetyl Coenzyme A carboxylase alpha around the phosphorylation site of Ser78 (PA(p-S)HK)
Target	Phospho-ACACA (Ser78)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	266 kDa
Uniprot ID	Q13085
Expiration Date	12 months from date of receipt.



Flow cytometry of HepG2(black line).Prim...