


## Product Datasheet

### Mouse SHH ELISA Kit (orb566111)

<b>Description</b>	Mouse SHH ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	SHH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q62226</a>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse SHH(Hedgehog Homolog, Sonic) ELISA...