


## Product Datasheet

### Mouse TBP2 ELISA Kit (orb566139)

<b>Description</b>	Mouse TBP2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	TBP2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<a href="#">Q6SJ95</a>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
TBP2(Thioredoxin  
Binding Protein  
2...