


## Product Datasheet

### Mouse TFR2 ELISA Kit (orb566149)

<b>Description</b>	Mouse TFR2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	TFR2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9JKX3</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
TFR2(Transferrin  
Receptor Protein  
...