


Product Datasheet

Mouse TPO-Ab ELISA Kit (orb566164)

Description	Mouse TPO-Ab (anti-Thyroid-Peroxidase Antibody) ELISA Kit
Reactivity	Mouse
Range	1.56-100IU/ml
Target	TPO-Ab
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938IU/ml
Expiration Date	Please enquire.



Mouse TPO-
Ab(anti-
Thyroid-
Peroxidase
Ant...