


Product Datasheet

Mouse alpha1-MG ELISA Kit (orb566184)

| | |
|------------------------|-------------------------------------------------|
| Description | Mouse alpha1-MG(alpha1-Microglobulin) ELISA Kit |
| Reactivity | Mouse |
| Range | 1.563-100ng/ml |
| Target | α 1-MG |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q07456 |
| Sensitivity | 0.938ng/ml |
| Expiration Date | Please enquire. |



Mouse α 1-MG(α 1-Microglobulin) ELISA Ki...