


## Product Datasheet

### Mouse alpha-CTx ELISA Kit (orb566185)

<b>Description</b>	Mouse alpha-CTx(Alpha Crosslaps) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	78.125-5000pg/ml
<b>Target</b>	$\alpha$ -CTx
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sensitivity</b>	46.875pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  $\alpha$ -CTx(Alpha Crosslaps) ELISA Kit