


## Product Datasheet

### Mouse CCL6 ELISA Kit (orb566202)

<b>Description</b>	Mouse CCL6 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	CCL6
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
CCL6(C-C  
motif  
chemokine 6)  
ELISA ...