


## Product Datasheet

### Mouse 4-HNE ELISA Kit (orb566243)

<b>Description</b>	Mouse 4-HNE ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	4-HNE
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse 4-HNE(4-Hydroxynonenal) ELISA Kit