


## Product Datasheet

### Mouse CXCR2 ELISA Kit (orb566266)

<b>Description</b>	Mouse CXCR2 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	CXCR2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<b>P35343</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
CXCR2(C-X-C  
chemokine  
receptor  
typ...