


Product Datasheet

Mouse GnRH ELISA Kit (orb566269)

Description	Mouse GnRH ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	GnRH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
GnRH(Gonadotropin
Releasing
Hormon...