


## Product Datasheet

### Mouse Calprotectin ELISA Kit (orb566273)

<b>Description</b>	Mouse Calprotectin ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1.25-80ng/ml
<b>Target</b>	Calprotectin
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, fecal extracts, Other liquid samples
<b>Uniprot ID</b>	<b>P27005, P31725</b>
<b>Sensitivity</b>	0.75ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
Calprotectin(Calprotectin)  
ELISA K...