

## Product Datasheet

### Mouse IL-17alpha ELISA Kit (orb566275)

**Description** Mouse IL-17alpha(Interleukin 17 alpha) ELISA Kit

**Reactivity** Mouse

**Range** 7.813-500pg/ml

**Target** IL-17 $\alpha$

**Storage** 2-8 °C for 6 months


**Note** For research use only

**Sample Types** Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples

**Uniprot ID** [Q62386](#)

**Sensitivity** 4.69pg/ml

**Expiration Date** Please enquire.



Mouse IL-17 $\alpha$ (Interleukin 17 alpha) ELIS...