


Product Datasheet

Mouse PDIA4 ELISA Kit (orb566277)

Description	Mouse PDIA4 ELISA Kit
Reactivity	Mouse
Range	0.313-20ng/ml
Target	PDIA4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P08003
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Mouse PDIA4(Protein disulfide-isomerase ...