


Product Datasheet

Mouse Ang1-7 ELISA Kit (orb566285)

Description	Mouse Ang1-7 ELISA Kit
Reactivity	Mouse
Range	7.8-500pg/ml
Target	Ang1-7
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P11859
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Mouse Ang1-7(Angiotensin 1-7) ELISA Kit