


Product Datasheet

Mouse IL36G ELISA Kit (orb566328)

Description	Mouse IL36G ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	IL36G
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q8R460
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
IL36G(Interleukin-
36 gamma) ELISA
...