


Product Datasheet

Mouse Dkk3 ELISA Kit (orb566343)

Description	Mouse Dkk3 ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Dkk3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QUN9
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse Dkk3
(Dickkopf-related
protein 3) ...