


Product Datasheet

Monkey ECF ELISA Kit (orb566374)

Description	Monkey ECF ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/ml
Target	ECF/CCL11
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	A0A0D9R208
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Monkey
ECF/CCL11(Eosinophil
Chemotactic ...