


Product Datasheet

Monkey IFN-beta ELISA Kit (orb566382)

Description	Monkey IFN-beta ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/ml
Target	IFN- β
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	O77812
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Monkey IFN- β (Interferon Beta) ELISA Kit