


Product Datasheet

Monkey IL-6 ELISA Kit (orb566394)

Description	Monkey IL-6 ELISA Kit
Reactivity	Monkey
Range	31.25-2000pg/ml
Target	IL-6
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P79341
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Monkey IL-6(Interleukin 6) ELISA Kit