


Product Datasheet

ACK1 (phospho-Tyr284) antibody (orb5664)

Description	ACK1 (phospho-Tyr284) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Ack1 (DH(p-Y)VM)
Target	Phospho-Ack1 (Tyr284)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	115 kDa
Uniprot ID	Q07912
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...