


Product Datasheet

Porcine ADP ELISA Kit (orb566457)

Description	Porcine ADP/Acrp30 (Adiponectin) ELISA Kit
Reactivity	Porcine
Range	62.5-4000ng/ml
Target	ADP/Acrp30
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q6PP07
Sensitivity	37.5ng/ml
Expiration Date	Please enquire.



Porcine
ADP/Acrp30(Adiponectin)
ELISA Ki...