


Product Datasheet

Porcine ANG ELISA Kit (orb566460)

Description	Porcine ANG ELISA Kit
Reactivity	Porcine
Range	0.313-20ng/ml
Target	ANG
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P31346
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Porcine
ANG(Angiogenin)
ELISA Kit