


Product Datasheet

Porcine ApoE ELISA Kit (orb566466)

Description	Porcine ApoE ELISA Kit
Reactivity	Porcine
Range	1.25-80ng/ml
Target	ApoE
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P18650
Sensitivity	0.75ng/ml
Expiration Date	Please enquire.



Porcine ApoE(Apolipoprotein E) ELISA Kit