


## Product Datasheet

### Porcine BALP ELISA Kit (orb566469)

<b>Description</b>	Porcine BALP ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	BALP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>A0A287BSC3</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
BALP(Bone  
Alkaline  
Phosphatase)  
...