


Product Datasheet

Porcine bFGF ELISA Kit (orb566472)

Description	Porcine bFGF ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/ml
Target	bFGF/FGF2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9TRD1
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Porcine
bFGF/FGF2(Basic
Fibroblast
Growt...