


Product Datasheet

Porcine cTn-I ELISA Kit (orb566485)

Description	Porcine cTn-I ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/ml
Target	cTn-I/TNNI3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	A5X5T5
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Porcine cTn-I/TNNI3(cardiac Troponin I) ...